



# MEDICALZ®

## TABLA DE TALLAS DE LAS PRENDAS

NOTA: Se recomienda encarecidamente medir a los pacientes antes de la cirugía. Las medidas son en pulgadas.

ZBRA' S/014, S/021, S/024,									MAMMARY STRAP S/013		PRESS-LIFT AC/002		
U.S. Sizes	30	32	34	36	38	40	42	44	Sizes	Thorax circ. (in)	Sizes	Neck Circ. (in) <sup>0a</sup>	Head Circ. (in) <sup>0a</sup>
Under Breast circ. (in) <sup>3</sup>	26¼ - 28¼	28¼ - 30¼	30¼ - 32¼	32¼ - 34¼	34½ - 36¼	36½ - 38¼	38½ - 40¼	40½ - 42	0 - 2X-Small	31½ - 33½	Infant	8½-9¼	16-18¼
Chest circumference (in) <sup>2</sup>									1 - X-Small	33½ - 35½	Toddler	10¼-11½	19¼-21¼
Cup A	30¼ - 31	32¼ - 33	34¼ - 35	36¼ - 37	38¼ - 39	40¼ - 41	42¼ - 43	44¼ - 45	2 - Small	35½ - 37½	Adult 1	11¼-13¼	22-24¼
Cup B	31 - 32	33 - 34	35 - 36	37 - 38	39 - 40	41 - 42	43 - 44	45 - 46	3 - Medium	37½ - 39½	Adult 2	14 - 16 ¼	25¼-26¼
Cup C	32 - 32¼	34 - 34¼	36 - 36¼	38 - 38¼	40 - 40¼	42 - 42¼	44 - 44¼	46 - 46¼	4 - Large	39½ - 41¼	Adult 3	16½-18½	27¼-28¼
Cup D	32¼ - 33½	34¼ - 35½	36¼ - 37½	38¼ - 39½	40¼ - 41½	42¼ - 43½	44¼ - 45½	46¼ - 47½	5 - X-Large	41¼ - 43¼			
Cup E	33½ - 34¼	35½ - 36¼	37½ - 38¼	39½ - 40¼	41½ - 42¼	43½ - 44¼	45½ - 46¼	47½ - 48¼					

LIPO-PANTY, COMBI & COMBI SUIT EC/001, EC/002, EC/002-H, EC/003, EC/004, EC/005, EC/006, EC/008, EC/009, EC/010, EC/016, EC/017, EC/034, EC/035		
Sizes	Hip circ. (in) <sup>3</sup>	Waist circ. (in) <sup>4</sup>
1 - 2X-Small	32½ - 34¼	23¼ - 25½
2 - X-Small	35 - 36½	25½ - 27½
3 - Small	37 - 39	27½ - 29½
4 - Medium	39½ - 41¼	29½ - 31½
5 - Large	41¼ - 46¼	31½ - 36¼
6 - X-Large	47¼ - 51¼	37 - 42¼
7 - 2X-Large	51½ - 55½	42½ - 47¼
8 - 3X-Large	56 - 59	47½ - 51¼
9 - 4X-Large	59½ - 62½	51½ - 55

BODY & BELT ELEGANCE COOLMAX EC/007, EC/011, EC/012, EC/033		
Sizes	Hip circ. (in) <sup>5</sup>	Waist circ. (in) <sup>4</sup>
1 - 2X-Small	32½ - 34¼	26¼ - 28¼
2 - X-Small	35 - 36½	28¼ - 30¼
3 - Small	37 - 39	31 - 33
4 - Medium	39½ - 41¼	33½ - 35¼
5 - Large	41¼ - 46¼	36 - 40
6 - X-Large	47¼ - 51¼	40½ - 44
7 - 2X-Large	51½ - 55½	44½ - 48
8 - 3X-Large	56 - 59	48½ - 52
9 - 4X-Large	59½ - 62½	52½ - 56

ABDOMINAL BINDER (Male or Female) S/022, S/22-E		
Sizes	Waist circ. (in) <sup>4</sup>	
1 - 2X-Small	25½ - 29½	
2 - X-Small	29½ - 33½	
3 - Small	33½ - 37½	
4 - Medium	37½ - 41¼	
5 - Large	41¼ - 45¼	
6 - X-Large	45¼ - 49¼	
7 - 2X-Large	49¼ - 53	
8 - 3X-Large	53 - 57	
9 - 4X-Large	57 - 61	
Height : 9 or 12.5 inches		

COMBI FOR MEN EC/037		
Sizes	Under Chest circ. (in) <sup>3</sup>	Hip circ. (in) <sup>5</sup>
1 - X-Small	37½ - 39½	35½ - 37½
2 - Small	39½ - 41¼	38¼ - 39¼
3 - Medium	41¼ - 43¼	40 - 42
4 - Large	43¼ - 45¼	42½ - 44½
5 - X-Large	45¼ - 47¼	44¼ - 46½
6 - 2X-Large	47¼ - 49¼	46¼ - 49

ROMEIO FOR MEN EC/038		
Sizes	Under Chest circ. (in) <sup>3</sup>	
0 - 2X-Small	35½ - 37½	
1 - X-Small	37½ - 39½	
2 - Small	39½ - 41¼	
3 - Medium	41¼ - 43¼	
4 - Large	43¼ - 45¼	
5 - X-Large	45¼ - 47¼	
6 - 2X-Large	47¼ - 49¼	

ABDOMINAL BINDER FOR MEN S/015		
Sizes	Hip circ. (in) <sup>5</sup>	
1 - X-Small	35½ - 37½	
2 - Small	38¼ - 39¼	
3 - Medium	40 - 42	
4 - Large	42½ - 44½	
5 - X-Large	44¼ - 46½	
6 - 2X-Large	46¼ - 49	

ARM SLEEVES & ADDTL MEASUREMENTS FOR COMBI EC/016, EC/042			
Sizes	Biceps circ. (in) <sup>7</sup>	Under Breast circ. (in) <sup>3</sup>	
1 - X-Small	9¼ - 11½	27½ - 31½	
2 - Small	11½ - 13	31½ - 35½	
3 - Medium	13 - 14½	35½ - 37½	
4 - Large	14½ - 16¼	37½ - 39½	
5 - X-Large	16¼ - 17¼	39½ - 41¼	
6 - 2X-Large	17¼ - 19¼	41¼ - 46	



Información y Ventas

Guanajuato No. 20 Interior A, Col. Roma C.P. 06700, CDMX  
T. (55) 55150590 y 57619378 / [cicamed@telmexmail.com](mailto:cicamed@telmexmail.com)  
[www.cicamed.com.mx/medicalz](http://www.cicamed.com.mx/medicalz)